

**matthew s yeh · bds · ms · msph**  
**specialist in orthodontics**

Bellingham  
(360) 676-9222

Stanwood  
(360) 629-4500

Lynden  
(360) 354-1349

**DR. YEH'S SECRET FORMULA**  
**FOR**  
**HOMEMADE APPLIANCE CLEANER**

This cleaner is far superior to commercial products, because it sterilizes your appliance *and* is more economical.

1. Brush appliance with toothbrush and water.
2. Soak appliance 20 minutes in the following solution:
  - 6 oz. water
  - 1 tsp. Chlorox
  - 2 tsp. Calgon (SPA BATH)
3. Rinse appliance with water and then mouthwash.